

## Standard Fixed Height Exam Couch - Navy Blue

**TASK**  
MEDICAL

The Task Medical Standard Fixed-Height Examination Couch is a heavy-duty, two-section couch with an easy hand-lift head section. Designed for durability and comfort, it provides long-lasting, maintenance-free service.

### Product Features

- Fixed height: 82 cm
- Heavy-duty round tubular steel legs for stability
- Easy hand-lift backrest with 10-stage angle adjustment
- Heavy-grade steel frame with a quality powder-coated finish
- Foam density: 36g/m<sup>2</sup>
- Hospital-grade upholstery: Fire, stain, mildew, oil, water, and wear-resistant PVC vinyl
- Also Available In: Grey & Black
- Australian-designed

CODE  
TKFHNB



### Product Specifications

- Total length: 184 cm
- Width: 71 cm
- Head section length: 72 cm
- Middle-leg section length: 108 cm
- Safe working load: 260 kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

